	ZIPPIE IRIS January 2018 EIZ5A								
			(()		EIZ5A				
	Date:						SUBMITTING FOR: Quote D	Order	
	Dealer Acct #:							01001	
	Dealer:						PO#:		
D	ealer Contact:					Sh	nip To:		
De	ealer Address:					Atte	ention:		
	Dealer City:		ST:ZIP:			Ade	dress:		
[Dealer Phone:	()	Fax: ()				dress:		
						Ship To	·	ZIP:	
	Confirm Via:	□ Fax				Ship To P	Phone: () Fax: ()		
	Mark For:								
		Thous	The HCPCS CODES provided are based				r interpretation of the code definitions. , verify the appropriate payer's coverage policy.		
					. For cover	-			
ŕ	Option #	HCPCS	Definition	Price	FRAME	Option #	HCPCS Definition	Price	
			Base Model	1.		Tre	moit (For use when Trendit Ontional forms is calestal)		
		М	laximum user weight is 225 lbs.		ᅵ┣	188TO45	ansit (For use when Transit Optional frame is selected) K0108 Transit Optional	\$250	
	Max		ght capacity for Z-Finity footrests is 165 lbs.				Tilt Release Options	+===	
	188WM48	E1233	Zippie IRIS Transit Included Includes True Fit Growth	\$2,815		188TT8	Foot Release N/A with Vent and Battery Tray.	NC	
	188WM47	E1233	Zippie IRIS Transit Optional Includes True Fit Growth	\$2,665	-		May limit folded width, folding frame defaults to right side □ Left □ Right (facing from rear)		
	188WM91	E1234	Zippie IRIS Folding Transit Included Includes True Fit Growth. N/A with Mono Back or Rocker Cover	\$2,940		188TT12	Dual Trigger Cable Release	NC	
		Zippie IRIS Folding Transit Optional			Pre-Set Tilt Angle (55° range) (Default in Bold, N/A w/ reclining back			g back)	
	188WM90	E1234	Includes True Fit Growth. N/A with Mono Back or Rocker Cover	\$2,790		188TA5	-5° to 50° N/A with 12.5" STFH	NC	
	Note: HCF	PCS Codes	s shown as submitted. Some sizes fit other H	CPC		188TA6	0° to 55°	STD	
		cc	odes i.e. 15-18" widths = E1161		, 0	188TA7	+5° to 60°	NC	
		Fra	me Struts (Rigid Frames Only)						
	188GY01		Fixed Struts Growth thru kits N/A w/ Folding Frame	STD		188TA8	Pre-Set Tilt Angle (40° range) (Default in Bold) -5° to 35° N/A with 12.5" STFH	NC	
			с С		- 	188TA9	0° to 40°	STD	
	188GY02		Width Adjusting System Built in growth of 2"	\$45		188TA10	+5° to 45°	NC	
			N/A w/ Mono Back or Folding		_				
			Frame Options						
	188FT46		Contracture Seat Frame N/A w/10" Frame Depth or Folding Frame	NC	_				
	188FS29		Battery Tray N/A w/ 10" & 11" frame widths, foot wheel lock or Folding Frame	\$295					
				2. SIZE	/ SEATII	NG			
		lidth (Cus	stom set-up on 17" & 18" folding frame thru B4				Custom Seating		
	188W1		10" - N/A w/ Mono Backrest	NC			stom seating to wheelchair integration, utilize universal h	ardware.	
	188W2		11" - N/A w/ Mono Backrest	NC NC		188QC01 188QC02	JAY SureFit		
	188W3 188W4		12" 13"	NC		188QC02 188QC03	JAY ConfigureFit		
	188W5		14"	STD					
	188W6		15"	NC			JAY Cushions (Optional)		
	188W7		16" - Folding Frame Grows to 18"	NC	To		cushions on-chair, please see the JAY price list. Cushio	n name or part	
	188W8		17" - Folding Frame Grows to 18"	NC	_	n	number must be filled out with desired width and depth.		
	188W9		18" - No Growth w/ Folding Frame	NC	1 -	400-1000	JAY Bundled Cushion		
			Seats			188JC20	Name / Part #	Price List	
	188ST35	E2231	Seat Pan-Growable Some limits may exist. See Page 9 for matrix of offering	\$200			Width: Depth:	1	
	188ST36	E2231	-	\$200	-				
	188ST54		Omit Seat Pan	STD	_				

2. SIZE / SEATING CONT.

Frame Depth									
				(Choose Select	ions fro	m A OR B)			
	Selection A:	Frame Dep	th Set up By Z	lippie			Selection B: Frame Depth Range		
	D	esired Seat	Depth			188SN01	Short Seat Frame 10"-14"	NC	
□ 10"	□ 11"	□ 12"	□ 13"	□ 14"		188SN02	Medium Seat Frame 12"-16"	NC	
□ 15"	□ 16"	□ 17"	□ 18"			188SN03	Long Seat Frame 14"-18"	NC	
□ 3" □ Flush	□ 2" □ 1" □ -1" □ -2"	Back Typ	Baci 1 3" fo	Type ward (+3')	OR B				
			3 1" fo 4	rward (+1') Flush		Selection B: Frame Depth			
This s	This section determines frame length					188FD1	10" rec. w/ Short Seat Frame N/A with Contracture Footrest System	NC	
					from	188FD2	11" rec. w/ Short Seat Frame	NC	
	Тор о	f Cushion to	Seat Rail		□ 5	188FD3	12" rec. w/ Short Seat Frame	NC	
□ 2 1/2" or 0	Greater			8		188FD4	13" rec. w/ Med. Seat Frame	NC	
Less that	n 2 1/2" (STD)				° □	188FD5	14" rec. w/ Med. Seat Frame	NC	
Below Se	Below Seat Rail				188FD6	15" rec. w/ Med. Seat Frame	NC		
			Γ	2 1/2" B		188FD7	16" rec. w/ Long Seat Frame	NC	
			E			188FD8	17" rec. w/ Long Seat Frame	NC	
				a france and a second		188FD9	18" rec. w/ Long Seat Frame	NC	

3. CASTERS / SEAT HEIGHT						
		Caster Wheel			Please cal	
188CW31		4" Micro Light Up - N/A with Folding	\$70		Rear Wh	
188CW73		4" x 1.25" Semi-Pneumatic -N/A w/ Folding	NC		Size	
188CW51		5" Micro Light Up	\$85			
188CW5		5" Low-Profile Polyurethane	STD		□ 12	
188CW7	E2219	5" x 1.5" Semi-Pneumatic	\$85		□ 12	
188CW33	E2214 E2215	6" Pneumatic	\$75			
188CW34		6" Polyurethane	NC		□ 12	
188CW35	E2219	6" x 1.5" Semi-Pneumatic	\$85			
188CW21	E2214 E2215	8" Pneumatic	\$85		Low Po	
188CW32		8" Polyurethane	NC			
188CW39	E2219	8" x 1.5" Semi-Pneumatic	\$85		□ 16	
188CW36	E2214 E2215	8" x 2" Pneumatic	\$145			
188CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert	\$225			

EAT HEIGHT								
Please call custo	mer servid	e for 13.2	5" on a 12'	' RW and 4	" Caster (Rigid	only)		
Rear Wheel	CASTER SIZE							
Size	□ 4"	□ 5"	□ 6"	8"x1"	□ 8"x1.5"	8"x2"		
	□ 13.0"	16.0"	16.0"	N/A	N/A	N/A		
□ 12"		16.5"	16.5"	N/A	16.5"	IN/A		
	N/A	17.0"	17.0"	17.0"	17.0"	17.0"		
		17.5"	N/A	17.5"	17.5"			
		15.0"	15.0"					
□ 12"	N/A	□ 15.5"	□ 15.5"	N/A	N/A	N/A		
		16.0"	16.0"					
Low Poly		16.5"	N/A	16.5"	16.5"	16.5"		
		13.5"	N/A	N/A	N/A			
		14.0"	14.5"			N/A		
□ 16"	N/A	□ 14.5" 	□ 17.0" 	□ 17.5" —	□ 17.5" —	_		
_		□ 17.5" —	□ 17.5" 	18.0"	□ 18.0" 	18.0"		
		18.0"	□ 18.0" □	18.5"	18.5"	_		
		18.5"	18.5"	N/A	□ 19.0"	□ 19.0"		
	N1/A	15.5"	□ 16.0"	N/A	N/A □ 16.5"	N/A		
□ 20"	N/A	□ 16.0" □ 16.5"	□ 16.5" □ 17.0"	□ 16.5" □ 17.0"	□ 16.5" □ 17.0"	□ 16.5" □ 17.0"		
		16.5"	□ 17.0" □ 16.5"	N/A	□ 17.0 □ 16.5"	N/A		
		□ 17.5"	□ 18.3 □ 17.0"	17.0"	□ 10.0"	17.0"		
□ 22"	N/A	□ 17.5 □ 18.0"	17.5	17.5	□ 17.5"	N/A		
		N/A	□ 18.0"	18.0"	□ 18.0"	□ 18.0"		
		N/A	17.0"	N/A		N/A		
		17.5"	17.5"	18.0"	18.0"	□ 18.0"		
□ 24"	N/A	18.0"	18.0"	18.5"	18.5"	N/A		
		18.5"	18.5"	19.0"	19.0"	19.0"		
N/A wit	h Foldin	g						

One Arm Drive
Only available with 6 spoke Mag wheels. One Arm Drive N/A with folding frame

188AL17	E0958	One Arm Drive Avail. w/ 20", 22", 24" Mag Rear Wheels	\$850
1880A2		Standard Handrim	NC
188HM1		Right	NC
188HM2		Left	NC
1880A1		Small Handrim	NC
188HM1		Right	NC
188HM2		Left	NC

Axle					
188AX26	One Arm Drive Axle Only avail. w/One Arm Drive Option	STD			
188AX1	Quick Release	NC			
188AX11	Quad Quick Release N/A with Lite Spokes	\$125			

	Rear Wheel Size						
12"-1	6" N/A w/handrii	ns. 12" w/attendant wheel lock N/A w/ +5 - 60° tilt angle.					
	188RW13	12" Mag - N/A with handrims	STD				
	188RW14	16" Mag - N/A with handrims	NC				
	188RW3	20" Mag	NC				
	188RW8	22" Mag	NC				
	188RW1	24" Mag	NC				
	188RW90	20" Lite Spoke	\$200				
	188RW91	22" Lite Spoke	\$200				
	188RW54	24" Lite Spoke	\$100				
	188RW103	20" Mag Wheel 5 Spoke	NC				
	188RW104	22" Mag Wheel 5 Spoke	NC				
	188RW102	24" Mag Wheel 5 Spoke	NC				

Rear Wheel Tire						
188RT1		Pneumatic	NC			
188RT2	E2211 E2213	Pneumatic Airless Insert	\$80			
188RT3		Full Polyurethane N/A w/ 12" & 16" rear wheels	NC			
188RT4		Low Polyurethane Only Available with 12" wheels	NC			

Rear Tire & Wheel Matrix

	Rear Wheel				
Tire	5 Spoke Mag	Spoke	Mag Wheel		
Pneumatic	20,22,24	24	12,16,20,22,24		
Pneumatic w/ Airless Insert	20,22,24	24	12,16,20,22,24		
Full Poly	20,22,24	24	20,22,24		
Low Poly	N/A	N/A	12		

12" = mm 16" = mm 20" = 451mm 22" = 501mm 24" = 540mm

	Handrim					
12" a	12" and 16" wheels not available with hand rims. See matrix for availability					
	188HR1	Aluminum Anodized	NC			
	188HR4	Plastic Coated	\$95			
	188HR7	Omit N/A with wheel/axle credit	NC			

Projections (Not Available on Plastic Coated)						
Only available with 24" Mag wheels & not available with spoke guards.						
188HR8	8 Oblique	\$325				

Rear Wheel & Handrim Matrix

	Rear Wheel				
Handrim	5 Spoke	Spoke	Mag		
	Mag	эроке	Wheel		
Aluminum Anodized	20,22,24	24	20,22,24		
Plastic Coat	20,22,24	24	20,22,24		
Projection 8 Oblique	N/A	N/A	24		

	5	. WHEEL LOO	CKS/HU	B LOCK		
Wheel Lock					Extension Handle	
Not available with Hub Lock				188EH1	E0961 6" Pair	\$75
188WL34	High Mount Push	NC				
40034/1.00	Foot Lock	NO		Only	Hub Lock	m drivo
188WL36		NC	- ⁻	,		\$100
188\//110	, ,	NC				\$100
	188WL36	Wheel Lock Not available with Hub Lock 188WL34 High Mount Push Foot Lock Foot Lock	Wheel Lock Not available with Hub Lock 188WL34 High Mount Push NC 188WL36 Foot Lock N/A w/ One Arm Drive, Folding frame defaults to left side, may limit folded width NC	Wheel Lock Not available with Hub Lock □ 188WL34 High Mount Push NC 188WL36 Foot Lock □ N/A w/ One Arm Drive, Folding frame defaults to left side, may limit folded width NC	Wheel Lock Not available with Hub Lock □ 188EH1 188WL34 High Mount Push NC 188WL36 Foot Lock N/A w/ One Arm Drive, Folding frame defaults to left side, may limit folded width NC 0 188LK01	Not available with Hub Lock 188WL34 High Mount Push NC 188WL36 Foot Lock N/A w/ One Arm Drive, Folding frame defaults to left side, may limit folded width NC 188LK01 12" Hub Lock

			6. FOOTRESTS	(See Appendix on Pa
	Hanger Mo	unted Fo	otrest (Footplates Mount to Hanger ⁻	Tubing)
-	Hanger	Type (Se	e Matrix on Page 9 for Lower Leg Length rai	nges)
	188H140		70° Swing In/Out Front Mount	NC
	188H142		80° Swing In/Out Front Mount	NC
	188H158		90° Swing In/Out Front Mount	NC
	188H9		70° Heavy Duty Front Mount Lift-Off	NC
	188H10		80° Heavy Duty Front Mount Lift-Off	NC
	188H11		90° Heavy Duty Front Mount Lift-Off	NC
	188H161	E0990	Elev.Fnt Mount - N/A 10-11" wide	\$295
	Exter	isions (m	ay interfere with some seat heights & caster	s)
	188E23		2" - not mounted	\$55
	188E24		4" - not mounted	\$55
			Footplate	
	188F9	K0040	Kids Angle Adjustable N/A with 17"-18" widths	\$185
	188F8	K0040	Adult Angle Adj. Available with 10"-18" widths	\$185
	188F62	K0040	Footboard Available with 10"-16" widths	\$150
			°	
	Ext. Tube	Mounted	Footrests (Footplates Mount to Ext.	Tube)
	-	Type (Se	e Matrix on Page 9 for Lower Leg Length ra	
	188H145		60° Swing In/Out	NC
	188H139		70° Swing In/Out	STD
	188H141		80° Swing In/Out - N/A with 8" Casters	STD
	188H159	K0053	Articulating - N/A 10"-13" widths	\$350
	188H160	E0990	Elevating N/A w/ platform footplates or < 14" width	\$295
	Extensi	ons (Se	e Matrix on Page 9 for Lower Leg Length ran	ges)
	188E1		Short	NC
	188E2		Medium	NC
	188E5		Long	NC
	188E16		High Mount	\$55
			Footplate	
	188F1		Composite N/A In 10", 11" or 12" width	STD
	188F3	K0040	Adult Angle Adjustable Available with 10"-18" widths	\$185
	188F4	K0040	Kids Angle Adjustable N/A with 17"-18" widths	\$185
	188F5		Aluminum Available with elevating 14"-18" wide	NC
	188F48	K0040	Platform Flip-Up - Ang Adj Available on 60° 70° and 80° ext. mount hangers only. N/A w/Folding Frame	\$190

	Hanger Mounts				
1" Offset N/A with -5-50° and -5-35° preset tilt angle and					
possible interference with 80° & 90° hangers					
188MZ02 Standard Hanger Mount STD					
188MZ01	1" Offset Hanger Mount	\$30			

age	9 for	matrix of Leg L	ength rang	ges)						
Γ			F	Footrest Hanger Options						
-		188IG2		Neoprene Covers Only available for ext. mount hangers	\$85					
		188HL1	E0951	Heel Loops N/A with Z-Finity 10-14 widths	\$30					
		188HL5	K0038	Leg Strap	\$35					
		188A54	K0038	Extra Single Leg Strap	\$30					
		188A110	E0952	Toe Loop Only avail. w/ adult ang adj. footplates, N/A with Z-Finity footrest	\$40					
		188H32		Omit Hanger	-\$15					
		188F7		Omit Footplate and Ext. Tube	-\$15					
Γ			Z-Fin	ity Hangers/Footrest System						
		Мах	kimum weig	ght capacity for Z-Finity footrests is 165 lbs.						
_		Sc	ome interfe	rence may exist in certain configurations.						
_		188H162		Z-Finity Hanger	STD					
		Z-Finity	Footrest	Mount (Must Chose see Illustration on Pag	,					
		188DI22		Direct Hanger Mt/Fixed Frame - High	NC					
		188DI23		Direct Hanger Mt/Fixed Frame - Low	NC					
_		188DI24		Extension Tube Mount	NC					
L		nity Footres	Pivot A	rm (Must Choose, use matrix on Page 9						
		188DI20		Short 2"	NC					
_		188DI21		Long 3 1/4"	NC					
				Z-Finity Footplate						
	_	Some	e interfere	ence may exist in certain configurations						
		188F320	K0040	Multi Angle Adjustable - Pair N/A with 8"-9" widths	\$185					
		188F321	K0040	Angle Adjustable Platform Footboard	\$150					
		188F6		Omit Footplate	NC					
		Tube Length (When Choosing Ext. Tube	Mount)							
_		188E62		Hemi Short	NC					
		188E1		Short	NC					
		188E2		Medium	NC					
		188E5		Long	NC					

JAY FOOT BOXES

Only available with Swing Away, Lift-off Hangers, Adult Angle Adj Footplates and ALR/ELR. Single Foot Box pairs N/A with 10"-15" frame widths. Full Foot boxes N/A with Heel Loop or Leg Strap or in odd number sizes, Tall Height N/A w/ 14"-15" frame widths and ALR/ELR. N/A with Z-Finity.

		Single Foot Box	
188JX01	E0954	Standard Height - Left	\$225
188JX02	E0954	Tall Height - Left	\$225
188JX03	E0954	Standard Height - Right	\$225
188JX04	E0954	Tall Height - Right	\$225
188JX05	E0954	Standard Height - Pair	\$450
188JX06	E0954	Tall Height - Pair	\$450
OR		Full Foot Box	
188JX07	E0954	Full Foot box Std Height	\$350
188JX08	E0954	Full Foot box Tall Height	\$350
188JX09	E0954	Full Foot box w/ Divider Std Ht	\$400
188JX10	E0954	Full Foot box w/ Divider Tall Ht	\$400

6. FOOTRESTS CONT. (See Appendix for matrix)

Contracture	Footrest	System	(CFS)
-------------	----------	--------	-------

	N/A w/ Neoprene Impact Guards, Kids Angle Adj.	Aluminum
188H26	Contracture Footrest System \$325	
188FP01	CFS - Center Mount NC	
188FP02	CFS - Outside Mount NC	
188E1	SHORT NC	C
188EP1	10.5" Extension Height NC	C
188EP2	11.0" Extension Height NC	
188EP3	11.5" Extension Height NC	
188EP4	12.0" Extension Height NC	
188EP5	12.5" Extension Height NC	
188EP6	13.0" Extension Height NC	


m Fo	otrests, Dual Po	est Armrest or Folding Frame.	
	188E2	REGULAR	NC
	188EP5	12.5" Extension Height	NC
	188EP6	13.0" Extension Height	NC
	188EP7	13.5" Extension Height	NC
	188EP8	14.0" Extension Height	NC
	188EP9	14.5" Extension Height	NC
	188EP10	15.0" Extension Height	NC

	7. BACKREST / SEATING							
DYNAMIC BACKREST					JAY Backs (Optional)			
		101 Dynamic Backrest N/A w/ Reclining or Folding Frame. N/A in some configurations.			Ordering	I JAY backs on-chair,	See the JAY Price List. Requires F	Part #
	188BT101		\$1,195			JAY Bundled Ba	ickrest	
					188JB02	Backrest Name	or Part #	Price List
						Width:	Height:	1

	STROLLER HAN	DLE BACK PO	OSTS	- BACK	RESTS (Cho	ose A, B,	C or D)
	A. Fixed Height						(
188BP9	Low (19" Ht)	NC					
188BP10	Medium (22" Ht)	NC					C. Mone
188BP11	Tall (25" Ht)	NC	Α			Must seled	
188BP13	Removable Extension	\$85					N/A w/
 188RG1	Folding Adjustable Rigidizer	\$175		Sele	ect One		
TOOKGT	N/A w/ Mono Back	φ175			188MO1		Mono Ba

	B Height Adjustable		
	Rigid Back w/ Angle Adjust Push Handle		
188B170	Std. (16"-22" Height.)	NC	
188PH8	Std. Angle Adj. Push Handle Lgth	NC	в
188PH9	Long Angle Adj. Push Handle Length	NC	P
188B171	Tall (19"-25" Height.)	NC	
188PH8	Std. Angle Adj. Push Handle Lgth	NC	1
188PH9	Long Angle Adj. Push Handle Length	NC	

OR



			C. Mono Back System		
		Must seled	ct JAY J3 when ordering Mono Back. N/A w/ Folding Frame.		
Sele	ct One				
	188MO1		Mono Backrest System	\$395	
	188MO2		Mono Backrest System with Dynamic Backrest	\$995	
AND	Select desir	ed Push	Handle		
	188PH8		Std. Angle Adj. Push Handle Length	NC	
	188PH9		Long Angle Adj. Push Handle Length	NC	

OR

		0	R		
		D Reclin	ning Backrest		1
 Only a	vailable v		It Range. N/A w/Folding Frai	ne	•
188BT21	E1225	Reclining	Backrest	\$850	
		188B172	Medium (21" Height)	NC	
		188B173	Tall (24" Height)	NC	
188BP64		Height adj	ij stroller handles [;] range 19"- 23" " widths w/Dual Trigger tilt	\$150	U
		188PH8	Standard push handle	NC	
		188PH9	Long push handle	NC	

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com Web: www.sunrisemedical.com ©1/2018 MK-100012 REV L Page 5

7. BACKREST / SEATING CONT.

Mono Backrest Options ONLY

		N/A on Fol	lding Frame, Vent/Battery Tray Base(LTV or	nly), Reclining E
		JAY J3	Back for Mono Backrest System	
	188JB04		JAY J3 Back Mono	See OF
	188JB05	JAY J3 via JAY Your Way Se		
JAY	J3 Backrest	Contour		
	188JP06	E2613	J3 PA Back - 2.2" of Contour Depth	\$535
	188JP08	E2615	J3 PL Back - 3" of Contour Depth	\$620
	188JP09	E2620	J3 PDL Back - 6" of Contour Depth N/A w/ 12" width	\$745
	188JP10	E2620	J3 PDC Back - 6" of Contour Depth N/A w/ 12" width	\$745
JAY	J3 Backrest	Height	•	
	188JR07		UT Standard - 16.5" shell height	NC
-		1		

	_	
188JR08	UT Medium - 18" shell height N/A w/ 12" width	NC
188JR09	UT Tall - 19.5" shell height N/A w/ 12" width	NC
188JR10	SH Standard - 21" shell height N/A w/ 12" width	NC
188JR11	SH Medium - 22.5" shell height N/A w/ 12" width	NC
188JR12	SH Tall - 24" shell height N/A w/ 12" width	NC


JAY J3 Backrest Width

Some interference may exist w/ armrest when ordering a backrest width greater than 4" wider than the base.

188JQ06	12" Backrest Width Available w/UT Standard 16.5" shell ht	NC
188JQ01	14" Backrest Width	NC
188JQ02	16" Backrest Width	NC
188JQ03	18" Backrest Width	NC

JAY J3 Backrest Color

188GZ01	Black	NC
188GZ02	Match Frame- Standard Co	olors NC
188GZ03	Match Frame- Custom Col	or \$50



	JA	Y J3 Ba	ck for Mono Backrest System Cont.	
Sele	ect desired He	adrest N	lount - optional	
	188MH10	E1028	Axys Mount & Bracket	\$230
	188MH15	E1028	Cobra Flip Back Mount and Bracket	\$280
	188MH16	E1028	Onyx Mount and Bracket	\$200
	188MH17		Omit	NC
Sele	ect desired He	adrest P	ad - optional	
	188HE65	E0955	Plush Head Support Pad 8"	\$230
	188HE66	E0955	Plush Head Support Pad 10"	\$230
	188HE78	E0955	Plush Head Support Pad 14"	\$230
	188HE87	E0955	HeadsUp Lrg: 5.3" Lng Arms	\$310
	188HE88	E0955	HeadsUp Med: 4.1" Sht Arms	\$310
	188HE89	E0955	HeadsUp Sml: 3.5" Sht Arms	\$310
	188HE83	E0955	Adjust-A-Plush Standard	\$260
	188HE82	E0955	Adjust-A-Plush Narrow	\$260
Sele	ect desired La	teral - op	tional	
	LAT 4X4 S/A	MOUNTS	S includes 1" OFFSET - PR	\$80
	188ZV14	E0956	LATERAL 4X4 W/ 1" OFFSET-PAIR	\$260
	188ZS07	E1028	LATERAL SWING AWAY MOUNT-PAIR	\$54
	LAT 4X5 S/A	MOUNTS	6 includes 1" OFFSET - PR	\$80
	188ZV15	E0956	LATERAL 4X5 W/ 1" OFFSET-PAIR	\$260
	188ZS07	E1028	LATERAL SWING AWAY MOUNT-PAIR	\$545
	LAT 4X6 S/A	MOUNTS	S includes 1" OFFSET - PR	\$80
	188ZV16	E0956	LATERAL 4X6 W/ 1" OFFSET-PAIR	\$260
	188ZS07	E1028	LATERAL SWING AWAY MOUNT-PAIR	\$54
	LAT 6X4 S/A	MOUNTS	6 includes 1" OFFSET - PR	\$805
	188ZV17	E0956	LATERAL 6X4 W/ 1" OFFSET-PAIR	\$260
	188ZS07	E1028	LATERAL SWING AWAY MOUNT-PAIR	\$54
	LAT 6X5 S/A	MOUNTS	S includes 1" OFFSET - PR	\$80
	188ZV18	E0956	LATERAL 6X5 W/ 1" OFFSET-PAIR	\$26
	188ZS07	E1028	LATERAL SWING AWAY MOUNT-PAIR	\$54
	LAT 6X6 S/A	MOUNTS	S includes 1" OFFSET - PR	\$80
	188ZV19	E0956	LATERAL 6X6 W/ 1" OFFSET-PAIR	\$26

Armrest - Height Adjustable

Dual Post not available with 0 - 13" deep or with 14 - 15" deep seat pan, Dual Post also not available with contracture Footrest. Some interference may exist in some configurations.

188AA12	E0973	T-Post Regular w/Foam	\$265
188AA13	E0973	T-Post Tall w/Foam	\$265
188AA3		Flip-Up	NC
188AA1	E0973	Adjustable Locking Flip-Up	\$255
188AA7	E0973	Length Adjust. Locking Flip-Up N/A with Contour Armpad	\$265
188AA62		Dual Post	NC
188AA63	E0973	Dual Post Height Adjustable	\$270
188AA43		Dual Post Flip Back	\$95
188AA48	E0973	Dual Post Flip Back Height Adj.	\$295
188AA11	E0973	Single-Post Height Adj. Std.	\$285
188AA29	E0973	Single-Post Height Adj. Low	\$285
188AA154	E0973	Height Adj. Cantilever Armrest Avail. w/ Mono Back only. N/A w/ 20", 22" and 24" rear wheels.	\$250

Armpad				
 Contour & Ergo only w/ Single Post Height Adj Foam only Length Adj.				
188AM23	Desk Classic	NC		
188AM21	Desk Contour	NC		
188AM24	Full Classic	NC		
188AM22	Full Contour	NC		
188AM9	Foam Short	NC		
188AM10	Foam Long	NC		
188AM27	Ergo armpad	NC		

Armrest F	Receiver	for Single Post Height Adjustable Arn	nrest
188AR17		Low	NC
188AR18		Std	NC

188JL01	E2209	Arm Trough - 14" long - Left	\$205
188JL04	E2209	Arm Trough - 14" long - Right	\$205
188JL05	E2209	Arm Trough - 14" long - Pair	\$410
188JL02	E2209	Arm Trough - 18" long - Left	\$205
188JL06	E2209	Arm Trough - 18" long - Right	\$205
188JL07	E2209	Arm Trough - 18" long - Pair	\$410
188JL03	E2209	Arm Trough - 22" long - Left	\$205
188JL08	E2209	Arm Trough - 22" long - Right	\$205
188JL09	E2209	Arm Trough - 22" long - Pair	\$410
	JA	Y Arm Trough - Hardware	
188JN01		Fixed - Each	NC
188JN03		Fixed - Pair	NC
188JN02	E1028	Swing Away - Rotational - Each	\$250
188JN04	E1028	Swing Away - Rotational - Pair	\$500

		Side Guard	
 Not		with height adjustable armrests, T Posts or over armrests. N/A with Mono Back.	
188SG33	K0108	Plastic Std	\$140
188SG29	K0108	Plastic Low	\$140

	Armrest Matr	ix for Single			
Armrest Size			Clas	sic Pad	
Armrest Size Std Low	Transfer Tube Length D	Receiver	H Range	Pad I	ength L
		Receiver H Range Pad Desk Std 10-3/4" to 15-3/4" Low Pad Desk Std 10-3/4" to 13-3/4" 9" to 13" Low 10-1/4" Std 9" to 13" Low 10-1/4" Receiver Contour Pad Receiver H Range Pad Desk Std 10" to 15" Low 10" to 15" Low 8-1/4" to 13-1/4" Std 10-1/2"	Full		
C+d	8"	Std	10-3/4" to 15-3/4"		
310	0	Low	8-3/4" to 13-3/4"	10 1/4"	14"
Low	5-1/4"	Std	9" to 13"	Desk 10-1/4" tour Pad Pad I Desk 10-1/2"	
Low	5-1/4	Low	7-1/4" to 11-1/4"		
			Contour Pad		
L		Receiver	H Danas	Pad length L	
			H Range	Desk	Full
		Std	10" to 15"	10 1/0"	14-1/4"
		Low	8-1/4" to 13-1/4"		
	Std	8-1/2" to 12-1/2"	10-1/2	14-1/4	
C	TO H	Low	6-3/4" to 10-3/4"	· · · · · · · · ·	
			Ergo Pad		
		Receiver	H Range	Pad l	ength L
		Std	10-1/2" to 15-1/2"		
	-D	Low	8-3/4" to 13-3/4"		2/4"
		Std	9" to 13"	11-3/4"	
		ora			

JAY Arm Troughs

N/A with Length Adj. Locking Flip Up and Adj Locking Flip Up armrests and Single Post Height Adj Armrests. Full Length Classic Armpads will automatically be selected when choosing only one side Arm Trough.

_				9. ACCE	ESSORIE	S	
		Back F	Pack (Kids Backpack included)				
		All Zippie c	hairs ship with one Kids Backpack at NC		Need	l 2" back post cl	learance to l
	188A1		Utility Backpack - Black	\$65	seat	ing; please note	
	188A2		Additional Zippie Backpack	\$60			with Mo
						1888A06	K0108
			Seat Pouch			1888A07	K0108
□ .	188A8		Black	\$50			
			Miscellaneous				JAY
	400400		1	\$45		188JY01	E0960
<u> </u>	188A83		Tool Kit	-φ 4 5		188JY02	E0960
		Р	ositioning Belt (w/clamps)			188JY03	E0960
	188A112	E0978	1" Dual Pull Side Squeeze Padded	\$85		188JY04	E0960
	188A117	E0978	Belt 2" Aircraft Buckle	\$110		188JY05	E0960
	188A118		Belt 2" Aircraft Buckle Padded	\$120			
	188A814	E0978	JAY 2 Pt 1" Push Btn - 4.5" Pad	\$120		188JZ01	
	188ZZ25	E0978	JAY 2 Pt 1.5" Push Btn - No Pad	\$40		188JZ02	
	188ZZ26	E0978	JAY 2 Pt 1.5" Push Btn - 7.5" Pad	\$86		188JZ03	
	188ZZ27	E0978	JAY 2 Pt 2" Push Btn - No Pad	\$40		188JZ04	
	188ZZ28	E0978	JAY 2 Pt 2" Push Btn - 9.5" Pad	\$86			
		Vertical N	Iount Vent Tray-Back cane Mount				
	188LF12	E1029	for LTV Vent (limit 15lbs) includes 1 Utility	\$495			J.
_			Hook Must remove for folding.			188JA01	E0960
	188LF13	E1029	Dovetail Vertical Frame Mt Hardware for Trilogy Vent (limit 15lbs) includes 1 Utility	\$495		188JA02	E0960
			Hook - Must remove for folding.	<i><i><i>ϕ</i></i></i>		188JA03	E0960
٦.	188LF7		Utility Hook	\$25		188JA04	E0960
		Up	per Extremity Support Tray				
l	•		ull armrest or desk armrest pads, add 1" to se			Not ava	ailable with a
	determine be	est tray size	. N/A with JAY Arm Trough. Must remove for t	olding.		188A23	K0065
ם [188A29	E0950	X-Sm (Suggested 10"-12" Width)	\$255		188A303	K0065
٦.	188A30	E0950	Sm (Suggested 13"-14" Width)	\$255			
Ξ.	188A31	E0950	Med (Suggested 15"-17" Width)	\$255			
	188A32	E0950	Lrg (Suggested 18"-20" Width)	\$255		188AT1	E0971
				10. FRA		DR	
		î	Standard Colors				
-	188C140		Aztec Gold	NC		188C116	
	188C2		Black - N/A as base with Zebra	NC		188C151	
-	188C141		Black Cherry	NC		188C152	
		1				188C150	
	188C10		Black Opal	NC			
	188C10 188C38		Blue Opal	NC		188C149	
	188C10 188C38 188C36		Blue Opal Candy Blue	NC NC		188C149	
	188C10 188C38 188C36 188C14		Blue Opal Candy Blue Candy Purple	NC NC NC			
	188C10 188C38 188C36 188C14 188C35		Blue Opal Candy Blue Candy Purple Candy Red	NC NC NC NC		188C43	
	188C10 188C38 188C36 188C14 188C35 188C142		Blue Opal Candy Blue Candy Purple Candy Red Electric Blue	NC NC NC NC NC			
	188C10 188C38 188C36 188C14 188C35 188C142 188C143		Blue Opal Candy Blue Candy Purple Candy Red Electric Blue Evergreen	NC NC NC NC NC NC		188C43 188C42	
	188C10 188C38 188C36 188C14 188C35 188C142 188C143 188C126		Blue Opal Candy Blue Candy Purple Candy Red Electric Blue Evergreen Glow	NC NC NC NC NC NC NC		188C43	
	188C10 188C38 188C36 188C14 188C35 188C142 188C143 188C126 188C144		Blue Opal Candy Blue Candy Purple Candy Red Electric Blue Evergreen Glow Green Apple	NC NC NC NC NC NC NC		188C43 188C42	
	188C10 188C38 188C36 188C14 188C35 188C142 188C143 188C126 188C144 188C147		Blue Opal Candy Blue Candy Purple Candy Red Electric Blue Evergreen Glow Green Apple Hot Sparkle Pink	NC NC NC NC NC NC NC NC		188C43 188C42	
	188C10 188C38 188C36 188C14 188C35 188C142 188C143 188C126 188C144 188C147 188C148		Blue Opal Candy Blue Candy Purple Candy Red Electric Blue Evergreen Glow Green Apple Hot Sparkle Pink Mauve Pink	NC NC NC NC NC NC NC NC NC NC		188C43 188C42	
	188C10 188C38 188C36 188C14 188C35 188C142 188C143 188C126 188C144 188C147 188C148 188C103		Blue Opal Candy Blue Candy Purple Candy Red Electric Blue Evergreen Glow Green Apple Hot Sparkle Pink Mauve Pink Rootbeer	NC NC NC NC NC NC NC NC NC NC NC		188C43 188C42	
	188C10 188C38 188C36 188C14 188C35 188C142 188C143 188C126 188C144 188C147 188C148 188C103 188C114		Blue Opal Candy Blue Candy Purple Candy Red Electric Blue Evergreen Glow Green Apple Hot Sparkle Pink Mauve Pink Rootbeer Sparkle Silver	NC NC NC NC NC NC NC NC NC NC NC NC		188C43 188C42	
	188C10 188C38 188C36 188C14 188C35 188C142 188C143 188C126 188C144 188C147 188C148 188C103		Blue Opal Candy Blue Candy Purple Candy Red Electric Blue Evergreen Glow Green Apple Hot Sparkle Pink Mauve Pink Rootbeer	NC NC NC NC NC NC NC NC NC NC NC		188C43 188C42	

Canopy (SunShade)

o mount canopy hardware, may require next height based on erference may occur with some headrest configurations. N/A Iono Back. Must remove for folding.

1888A06	K0108	Canopy (Sunshade) 10"-14" Width	\$275
1888A07	K0108	Canopy (Sunshade) 15"-20" Width	\$275

JAY Anterior Thoracic Support									
Style									
188JY01	188JY01 E0960 Classic Dynamic								
188JY02	E0960	Classic Structure	\$187						
188JY03	E0960	Contour Dynamic	\$187						
188JY04	E0960	Contour Structure	\$187						
188JY05	E0960	Center Opening with Buckle	\$187						
		Size							
188JZ01		Small	NC						
188JZ02		Medium	NC						
188JZ03	NC								
188JZ04		X-Large	NC						

JAY Padded Trunk Strap									
	188JA01	E0960	\$187						
	188JA02	E0960	Padded Trunk Strap - Medium	\$187					
	188JA03	E0960	Padded Trunk Strap - Large	\$187					
	188JA04	E0960	E0960 Padded Trunk Strap - XL						
			Spoke Guard						
Not available with 12" & 16" rear wheels & Projections Handrims									
	188A23	K0065	Black	\$115					
	188A303	K0065	0065 Clear Zippie						

	Anti-Tip									
		188AT1	E0971	Rear Anti-Tip	\$140					
E (COLOR									
	Standard Matte Colors									
		188C116		Matte Black	NC					
		188C151		Matte Black Cherry	NC					
		188C152		Matte Purple	NC					
		188C150		Matte Electric Blue	NC					
		188C149		Matte Evergreen	NC					

Custom Colors								
188C43 Mossy Oak Camouflage \$29								
188C42	Stars and Stripes	\$295						
188C108	Zebra Choose base color from Standard colors: White is Standard.	\$295						



Seat Pan Size

10" - 16"

14" - 18"

Growing Seat Pan Matrix - Rigid and Folding Frame Options

Rigid Frame - Frame Depth Growth Range Matrix

000

000 0

0 0 Å.

DIM A

Ŧ

I De

000

000

0



Footrest Adaptor

Short

Long



E	Direct Hang	jer Mounte	d	E:	xtensio
	Hanger			Footrest	Extens
	Mounting	Min	Max	Adaptor	Tube
	Position			Auaptor	Leng
	High	-0.9"	3.4"		Hemi S
	Mount	-0.9	3.4	Short	Shor
	Low	2 6"	6.7"	Short	Mediu
	Mount	2.0	0.7		Long
	High	-2.1"	4.6"		Hemi S
	Mount	-2.1	4.0	Long	Shor
	Low	1.3"	7 9"	Long	Mediu
	Mount	1.3	1.9		Long

		Seat Pan to T		
Hanger	Extension	Min	Max	Min Allowable FSH
	No Ext.	3 1/2"	9"	11 1/2"
70° Front Mount	2"	3 1/2"	11 1/2"	13 1/2"
	4"	3 1/2"	13 1/2"	15 1/2"
	No Ext.	3 1/2"	9 1/2"	11 1/2"
80° Front Mount	2"	3 1/2"	12"	14"
	4"	3 1/2"	14"	16"
90° Front Mount	No Ext.	3"	9"	11 1/2"
ELR Front Mount	No Ext.	3"	7"	10 1/2"
	Short	9"	12 1/2"	12 1/2"
60° Swing In/Out	Medium	9"	14 1/2"	12 1/2"
	Long	11 1/2"	16"	14 1/2"
	Short	10 1/2"	15 1/2"	12 1/2"
70° Swing In/Out	Medium	10 1/2"	17 1/2"	13"
	Long	13 1/2"	19"	15"
	Short	11"	14 1/2"	13"
80° Swing In/Out	Medium	11"	16 1/2"	13"
	Long	13"	18"	15"
	Short	14"	17 1/2"	12 1/2"
ALR	Medium	15"	19 1/2"	13 1/2"
	Long	16 1/2"	19"	14 1/2"
	Short	13"	17"	12 1/2"
ELR Ext Mount	Medium	13"	18 1/2"	13 1/2"
	Long	14 1/2"	20"	14 1/2"





E	xtension Tu	ube Mounte	ed		
Footrest Adaptor	Extension Tube Length	Min	Max		
	Hemi Short	3.9"	12.4"		
Short	Short	4.9"	13.4"		
Short	Medium	6.6"	15.4"		
	Long	8.6"	17.4"		
	Hemi Short	2.6"	13.6"		
Long	Short	4.1"	14.6"		
Long	Medium	5.6"	16.6"		
	Long	7.6"	18.6"		

Folding Frame - Frame Depth Growth Range Matrix Seat Depth Seat Depth

Dim A Range

11-1/4" - 15-1/4"

15" - 19'

Seat Depth					500		Put					L Kange		
Pan Size	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	Min	Max	
Short	х	х	х	х	х							10"	14"	
Medium					х	х						11.75"	16.25"	
Long							х	х				13.75"	18.25"	
X-Long									х	х	x	15.75"	20.25"	

Depth

0 0 0 Ş 0 0 0 0

Additional Options Through Built-4-Me

To inquire about any unique customization for Zippie products, contact: Phone: (800) 333-4000, Fax: (800) 333-9011

built<mark>4</mark>me

Commonly Requested Options Include:

- √ **Reverse Configuration**
- Angle Adjustable Handle for Reclining Backrest
- ⁄ **Dual Mount Wheel Locks - Attendant Foot Lock**
- 1 Dynamic Back on Frame Depth up to 20"
- **O2 Holder with Transit** 1
- IV Pole with Transit 1
- √ Solid Seat w/ Quick Release Hardware for folding frame



Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com Web: www.sunrisemedical.com ©1/2018 MK-100012 REV L Page 9